



GynTect®

Epigenetic markers for cervical cancer diagnostics

Info for labs



- CE-IVD approved diagnostic test* identifying patients with malignant changes of the cervix
- Performance possible in one working day
- Performed on cobas® z 480 Analyzer
- Evaluation using common calculation software


* The product is currently only available in the following countries: Canada, Italy, Poland, Portugal and Turkey.



GynTect® – reliable and fast diagnostics

If the cervical cancer screening result is **abnormal** with the **Pap test** and/or **positive** with the **HPV test**, patients suddenly find themselves in an exceptional situation. Even though both tests indicate a possible cancer, **in many cases** there is **no malignant disease** and the positive test result was a false alarm. Further examinations are necessary for reliable clarification, such as a colposcopy with biopsy, if required. If abnormalities are detected, the presumably affected tissue is often removed.

GynTect® is a fast and non-invasive test for clarification of abnormalities in cervical cancer screening. Just one further smear allows for a reliable result to be obtained within a few days.



An existing infection with **HPV** may lead to **genetic instability of the infected cells** and eventually cervical cancer. In the course of **carcinogenesis, changes (methylations) occur in the DNA.**

GynTect® recognises six areas of the human genome, which only exist methylated during the development of cancer cells. **GynTect®** thus identifies patients with malignant cervical cell changes.

Decision-making based on reliable results

With a **negative GynTect®** result, a **cancer diagnosis** could be **excluded** at the time of testing. If there was an abnormal Pap test result or HPV infection present prior to the test, it is recommended to observe them further.

If there is a **positive GynTect®** result, a **malignant precursor or even cancer** is **very likely**. Further steps such as diagnostics assisted by colposcopy and surgical therapies are recommended.

Based on available study data, **GynTect®** provides a clear indication of malignancy status in patients with abnormal Pap smear: In all previous studies, **GynTect®** was able to detect all cases of cervical cancer (sensitivity = 100%).

GynTect® is rarely positive in patients with inconspicuous cytological findings (specificity = 96.6%). Cancer develops via the histopathologically defined dysplasias CIN1, CIN2 and CIN3. **GynTect®** detection rates for these dysplasias increase continuously. This indicates a prognostic value of the **GynTect®** cancer markers.



Negative
GynTect® result

Cervical cancer is very unlikely
at the time of testing.

If there is a dysplasia, it is very unlikely
to be malignant.




Positive
GynTect® result

Cervical cancer or a malignant precursor
is very likely.

Therapeutic action is recommended!

Study data GynTect®



Detection rate **GynTect®** (dark purple bar) depending on the clinical status of the patient

Confidence interval for confidence level = 95%:

NILM: 2.22–4.99%

CIN1: 7.96–34.16%

CIN2: 19.92%–46.32%

CIN3: 57.5–70.11%

Cancer: 88.06–100%



Assay principle and workflow

The **GynTect**[®] assay principle is based on the detection of **DNA methylation** in human gene regions that occurs **specifically during carcinogenesis**. In the process of DNA methylation, methyl groups are added to the DNA. These are always cytosines located next to guanines (**CpG dinucleotides**).

The analysis of a patient sample comprises two steps:

1. First, methylation is fixed by **bisulphite treatment**.
2. Subsequently, specific regions of the genome are analysed by **PCR** and an **evaluation** is carried out using common spreadsheet software.



Only originally methylated DNA regions are amplified in the PCR. Therefore, this procedure is also called **methylation-specific PCR (MSP)**.

The **GynTect**[®] assay includes **several internal controls** to ensure a highly reliable and robust workflow. Moreover, positive and negative controls are included.