




# MERGITE!

## Automated washing of IIFT slides



- ✓ **Fast, parallel washing** of 50 substrate fields
- ✓ **Washing each field separately**
- ✓ **Standardisation ensures high reproducibility of results**
- ✓ **No beakers, cuvettes, and paper tissue**
- ✓ **Intuitive user interface** on the integrated touch screen
- ✓ **Excellent service** from EUROIMMUN

Versions .....	MERGITE! 10 for 5 slides with 10 fields each, MERGITE! 50 for 1 slide with 50 fields
Power supply .....	100–240 V, 100 VA , 47/63 Hz
Width x depth x height .....	250 x 381 x 207 mm
Weight .....	MERGITE! 10: 7.5 kg, MERGITE! 50: 6.8 kg