





HPV microarray (IVD) – complete typing of human papilloma viruses



- Detection and typing of all 30 relevant anogenital HPV subtypes in one reaction
- Direct detection of the viral oncogenes E6/E7 provides highest possible sensitivity
- Significant results even in very early stages of the infection
- Distinction between high-risk and low-risk types of HPV
- Reliable identification of multiple infections

EUROArray HPV

(order no. MN 2540-####)



####: For detailed information about available test kit formats see our product catalogue or visit www.euroimmun.com.

- Well-established EUROArray technology
- Simple test performance – no in-depth molecular biology knowledge required!
- Ready-to-use PCR components, integrated controls
- Fully automated and standardised evaluation, interpretation and archiving of results (EUROArrayScan system)