




First multiparameter test for Mediterranean-specific inhalation allergens

EUROLINE DPA-Dx* Pollen Southern Europe 1

Component-based diagnostics for optimal therapy decisions



Differential diagnostics in a new dimension

- Clinically approved multiparameter test recommended by leading experts in the field
- Identification of the allergy-causing proteins, clinically relevant in Southern Europe (e.g. nCup a 1, rOle e 1, rPar j 2, rBet v 1, rPhl p 1, rPhl p 5)
- Differentiation between primary sensitisations and cross reactions
- Targeted selection of patients suitable for a specific immunotherapy

*Defined partial allergen diagnostics