




# EUROArray ArboCDZ\*

Direct detection and typing of chikungunya, dengue and Zika viruses in one reaction



## One PCR – seven parameters

Chikungunya virus

Dengue virus type 1  
Dengue virus type 2  
Dengue virus type 3  
Dengue virus type 4

Zika virus Asian lineage  
Zika virus African lineage



Reliable differentiation between seven different arboviruses (chikungunya, dengue, Zika)



Simple test procedure:

- Ready-to-use reagents
- Integrated controls
- Fully automated and standardised evaluation



Sample material:

- Serum/plasma
- Urine

\* For research use only, not for in vitro diagnostic use