


euroimmun 

From Revvity

Test systems for infection diagnostics.



Direct detection of infectious agents



Pathogen	Antigen ELISA	PCR	Microarray
<i>Aspergillus</i>	●		
<i>Chlamydia trachomatis</i>			● ³
Dengue-Viren-NS1	●		
Dermatophytes (n = 23) ¹ and yeast/moulds (n = 6) ²			●
<i>Haemophilus ducreyi</i>			● ³
Herpes simplex virus (HSV-1/2)			● ³
Human papillomavirus (HPV)			●
<i>Mycoplasma genitalium</i>			● ³
<i>Mycoplasma hominis</i>			● ³
<i>Neisseria gonorrhoeae</i>			● ³
<i>Toxoplasma gondii</i>		●	
<i>Treponema pallidum</i>			● ³
<i>Trichomonas vaginalis</i>			● ³
<i>Ureaplasma parvum</i>			● ³
<i>Ureaplasma urealyticum</i>			● ³

¹anthropophilic: *T. tonsurans*, *T. interdigitale*, *T. schoenleinii*, *T. concentricum*, *T. rubrum*, *T. violaceum*, *E. floccosum*, *M. ferrugineum*, *M. audouinii*; zoophilic: *T. equinum*, *T. mentagrophytes*, *T. simii*, *T. quinckeanum*, *T. erinacei*, *T. bullosum*, *T. benhamiae*, *T. verrucosum*, *T. eriotrephon*, *M. canis*, *N. persicolor*; geophilic: *N. fulva*, *N. gypsea*, *N. incurvata*;
²*C. albicans*, *C. guilliermondii*, *C. parapsilosis*, *F. oxysporum*, *F. solani*, *Sc. brevicauli*;
³included in the EUROArray STI (multiplex test)

Detection of antibodies against infectious agents



Pathogen	ELISA	Blot	IFA	ChLIA
Adenovirus	●			
<i>Bartonella</i>			●	
<i>Bordetella</i> FHA	●	●		
<i>Bordetella pertussis</i> toxin	●	●		
<i>Borrelia</i>	● ⁴	● ⁴	●	● ^(*)
<i>Brucella abortus</i>	●			
<i>Campylobacter jejuni</i>	●			
Chikungunya virus	●		●	
<i>Chlamydia pneumoniae</i>	●		●	
<i>Chlamydia psittaci</i>			●	
<i>Chlamydia trachomatis</i>	●		●	
<i>Coxiella burnetii</i> , phase 1	●			
<i>Coxiella burnetii</i> , phase 2	●			

* manufactured by Technogenetics. Please contact your local Euroimmun representative for more information.
^(*)manufactured by Euroimmun or Technogenetics;
⁴also CSF diagnostics;

Pathogen	ELISA	Blot	IFA	ChLIA
Coxsackievirus			●	
Crimean Congo fever virus	●			
Cytomegalovirus (CMV)	● ^{4,5}	●		● [*]
Dengue virus	●		●	
Diphtheria toxoid	●			
<i>Echinococcus</i>	●	●		
Enterovirus	●			
Epstein-Barr virus (EBV)	● ^{4,5}	●	● ⁵	● ^(*)
Hantavirus	●		●	
<i>Helicobacter pylori</i>	●	●		
Hepatitis E virus (HEV)	●	●		
Herpes simplex virus, type 1 (HSV-1)	●	●		● ^(*)
Herpes simplex virus, type 2 (HSV-2)	●	●		● ^(*)
Herpes simplex virus (HSV-1/2)	● ⁴	●		● [*]
Human herpesvirus 6 (HHV-6)			●	
Influenza virus, type A	●			
Influenza virus, type B	●			
Japanese encephalitis virus	●		●	
<i>Legionella pneumophila</i>	●		●	
<i>Leishmania</i>	●		●	
<i>Listeria monocytogenes</i>			●	
Mayaro virus	●			
Measles virus	● ^{4,5}			● [*]
MERS coronavirus (MERS-CoV)	●			
Mumps virus	● ⁴			● [*]
<i>Mycoplasma pneumoniae</i>	●			
Parainfluenza virus, types 1–4	●			
Parvovirus B19	●			● [*]
<i>Plasmodium</i>	●			
Respiratory syncytial virus (RSV)	●			
Rubella virus	● ^{4,5}			● [*]
Sandfly fever virus	● ⁶		●	
SARS coronavirus 2 (SARS-CoV-2)	●			●
<i>Schistosoma</i>	●		●	
<i>Strongyloides</i>	●			
TBE virus	● ⁴		●	
Tetanus toxoid	●			● [*]
<i>Toxocara canis</i>	●	●		
<i>Toxoplasma gondii</i>	● ^{4,5}		●	● [*]
<i>Treponema pallidum</i>	● ⁴	●	●	
<i>Trichinella</i>	●			

* manufactured by Technogenetics. Please contact your local Euroimmun representative for more information.
^(*)manufactured by Euroimmun or Technogenetics;
⁴also CSF diagnostics;
⁵also avidity determination;
⁶only Toscana virus available

Pathogen	ELISA	Blot	IFA	ChLIA
<i>Trypanosoma cruzi</i>	●			
Usutu virus	●			
Varicella zoster virus (VZV)	● ^{4,5}		●	● [*]
West Nile virus	● ⁵		●	
Yellow fever virus			●	
<i>Yersinia enterocolitica</i>	●	●		
Zika virus	●		●	

* manufactured by Technogenetics. Please contact your local Euroimmun representative for more information.

(*) manufactured by Euroimmun or Technogenetics;

⁴also CSF diagnostics;

⁵also avidity determination

Detection of other markers



Detection of	Antigen ELISA
CXCL13	●

Find out more about infection diagnostics at
www.euroimmun.com
 or contact us directly at
www.euroimmun.com/contact

