




# EUROArray STI (IVD)

Detection of up to 11 relevant sexually transmitted pathogens in one reaction

BACTERIAL PATHOGENS
Chlamydia trachomatis (CT)
Haemophilus ducreyi (HD)
Mycoplasma genitalium (MG)
Mycoplasma hominis (MH)
Neisseria gonorrhoeae (NG)
Treponema pallidum (TP)
Ureaplasma parvum (UP)
Ureaplasma urealyticum (UU)

PROTOZOON
Trichomonas vaginalis (TV)



VIRAL PATHOGENS
Herpes simplex virus 1 (HSV-1)
Herpes simplex virus 2 (HSV-2)



- **Reliable detection of 11 sexually transmitted pathogens**
- **Detection of multiple infections**
- **Suitable for pathogens that cannot be cultivated**
- **Possible prevention of infertility and abortion**

## EUROArray STI

(order no. MN 2830-####)



####: For detailed information about available test kit formats see our product catalogue or visit [www.euroimmun.com](http://www.euroimmun.com)

- **Comprehensive multiplex testing**
- **Fully automated evaluation**
- **Integrated controls for all work steps**
- **Simple procedure – just one PCR, ready-to-use reagents**
- **User-friendly EUROArray technology**