



# Immunoblot PreQ

Pre-equipped (PreQ) incubation trays for fully automated processing of EUROLINE blot strips



- ✓ **Lowest manual work load** owing to pre-equipment of the incubation channels
- ✓ **Reliable fixation** of the blot strips at the bottom of the channels
- ✓ **User-friendly click system** allows flexible combination of blot trays
- ✓ **Fully automated incubation and evaluation** by EUROBlotOne and EUROLineScan Software

## Currently available PreQ immunoblot products for autoimmune diagnostics

Order no.	Product
DL 1530-5001-4 G	EUROLINE Autoimmune Inflammatory Myopathies
DL 1532-5001 G	EUROLINE Systemic Sclerosis Profile
DL 1590-5001-3 G	EUROLINE ANA Profile 3
DL 1590-5001-8 G	EUROLINE ANA Profile 1
DL 1590-5001-30 G	EUROLINE ANA Profile 3 plus DFS70
DL 1111-5001-7 G	EUROLINE Paraneoplastic Neurological Syndromes 12 Ag
DL 1300-5001-4 G	EUROLINE Autoimmune Liver Diseases
DL 0160-5001 G	EUROLINE Validation
ZG 0150-0101-146	Positioning aid for EUROBlotOne (fixation of incubation channels)